

糞便檢查採集程序

Stool collection procedure

步驟 1(Step 1)

採集新鮮糞便，採檢容器為 S-Y 糞便反應瓶(S-Y feces test bottle)。
Please collect fresh stool into S-Y feces test bottle.

步驟 2(Step 2)

採集前請其參閱包裝袋上之說明。Please read the instruction before collection.

步驟3(Step 3)

請將糞便解於衛生紙上，腹瀉時可先以紙杯接取部分糞便。
Please collect stool from toilet papers or from paper cup in case of diarrhea.

步驟4(Step 4)

打開S-Y糞便瓶瓶蓋，以採便棒挖取花生米大小之糞便，放入瓶內後旋緊瓶蓋。Pick up stool as big as a peanut in the test bottle, then tighten the cover.

步驟5(Step 5)

若為液狀便，則倒入 1-2ml 即可。Pour 1-2 ml in case of diarrhea.

步驟 6(Step 6)

檢體採集後，採檢容器須清楚標示姓名與病歷號。After collect, please write your name and chart number on the bottle.

步驟 7(Step 7)

將檢體裝於檢體袋前袋。
Put the bottle into the front package.

步驟8(Step 8)

將檢體袋盡速送交二樓檢驗科。
Send the package to OPD Laboratory on 2F, building,as soon as possible.

※注意事項 Attention※

- 1.為避免動物血及其它東西干擾，故檢查前 2 天禁食含動物血食物如豬肝、豬血...等動物內臟或紅肉類.魚肉等白肉可食用少許。
- 2.Aspirin 因會引起腸胃出血的副作用引起偽陽性。
- 3.香蕉、蘿蔔類,因內含類過氧化酵素會引起偽陽性。
- 4.大量食用 Vitamin C 藥片則會有抑制類過氧化酵素會引起偽陰性。